


Barbed Sink Head




Features:

- Sink the hose. Commonly used to connect the external tube of the peristaltic pump to prevent the tubing from twisting, floating and the liquid splashing
- Plug-in type, simple & convenient to use. Anti-slip and catch-edge design for no leakage
- SUS304 stainless steel material, safe and hygienic, featured with high strength, excellent heat resistance, anti-corrosion and no magnetic

★ Application Note



The bottom of the barbed sink head is designed with a **120° groove**, which enables the liquid to flow out stably in time while maintaining its own stability.



Unit: mm

Model No.	Tubing ID	Dimension (mm)						Material
Barbed Sink Head								
CT-016	1.0~1.6 mm	A=Φ0.99	B=Φ1.60	C=Φ2.80	D=Φ15	E=1	F=10.7	SUS304
CT-024	1.6~2.4 mm	A=Φ1.58	B=Φ2.40	C=Φ3.40	D=Φ25	E=2	F=14	SUS304
CT-032	2.4~3.2 mm	A=Φ1.97	B=Φ3.20	C=Φ4.20	D=Φ25.7	E=2.5	F=15.2	SUS304
CT-040	3.2~4.0 mm	A=Φ2.96	B=Φ4.00	C=Φ6.00	D=Φ27.3	E=3	F=20	SUS304
CT-048	4.0~4.8 mm	A=Φ3.45	B=Φ4.80	C=Φ6.80	D=Φ35.6	E=3.5	F=23	SUS304
CT-064	4.8~6.4 mm	A=Φ4.93	B=Φ6.40	C=Φ8.40	D=Φ40.7	E=4	F=25.8	SUS304
CT-087	6.4~8.7 mm	A=Φ6.40	B=Φ8.70	C=Φ10.70	D=Φ45.2	E=5	F=27	SUS304